

UFC1 RABBIT PAB

Cat.#: S218320

Product Name: Anti-UFC1 Rabbit Polyclonal Antibody

Synonyms: NEDSG; HSPC155

UNIPROT ID: Q9Y3C8 (Gene Accession - BC005187)

Background: UFC1 is an E2-like conjugating enzyme for ubiquitin-fold modifier-1 (UFMI; MIM 610553) (Komatsu et al., 2004 [PubMed 15071506]).[supplied by OMIM, Mar 2008]

Immunogen: Fusion protein of human UFC1

Applications: ELISA, IHC

Recommended Dilutions: IHC: 200-400; ELISA: 5000-10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human, Mouse, Rat


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin embedded Human liver cancer tissue using 218320(UFC1 Antibody) at a dilution of 1/200(Secreted).



In comparison with the IHC on the left, the same paraffin-embedded Human liver cancer tissue is first treated with the fusion protein and then with 218320(Anti-UFC1 Antibody) at dilution 1/200.