

UPK3B RABBIT PAB

Cat.#: S210023

Product Name: Anti-UPK3B Rabbit Polyclonal Antibody

Synonyms: FLJ32198, MGC10902, p35, UPIIIb

UNIPROT ID: Q9BT76 (Gene Accession - BC004304)

Background: UPK3B is a minor component of the apical plaques of mammalian urothelium that binds and dimerizes with uroplakin-1b (UPK1B, MIM 602380), one of the major conserved urothelium membrane proteins. The other major conserved integral membrane proteins of urothelial plaques are UPK1A (MIM 611557), UPK2 (MIM 611558), and UPK3A (MIM 611559)

Immunogen: Fusion protein of human UPK3B

Applications: ELISA, WB, IHC

Recommended Dilutions: IHC: 50-200;WB: 1000-4000;ELISA: 500-10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Signal Transduction


Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human brain tissue using 210023(UPK3B Antibody) at a dilution of 1/25(Cytoplasm or Nucleus).



Immunohistochemistry analysis of paraffin-embedded Human tonsil tissue using 210023(Anti-UPK3B Antibody) at a dilution of 1/25.



Gel: 10%SDS-PAGE, Lysate: 40 µg;
Lane: 293T cells;
Primary antibody: 210023(UPK3B Antibody) at dilution 1/1000;
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution;
Exposure time: 1 minute