

VILLIN RABBIT MAB

Cat.#: N263069

Product Name: Anti-Villin Rabbit Monoclonal Antibody

Synonyms: VIL; VIL1; Villin1

UNIPROT ID: P09327

Background: Ca(2+)-regulated actin-binding protein.

Immunogen: A synthetic peptide of human Villin

Applications: WB,IHC-P,IP

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R07-4A2

MW: Calculated MW: 93 kDa; Observed MW: 93 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Villin in MCF-7 lysates using Villin antibody.