

ZFP64 RABBIT PAB

Cat.#: S218044

Product Name: Anti-ZFP64 Rabbit Polyclonal Antibody

Synonyms: ZNF338

UNIPROT ID: Q9NTW7 (Gene Accession - BC021087)

Background: Predicted to enable DNA binding activity and metal ion binding activity. Predicted to be involved in positive regulation of cytokine production and positive regulation of transcription by RNA polymerase II. Predicted to act upstream of or within positive regulation of mRNA splicing, via spliceosome. Predicted to be active in nucleus.

Immunogen: Fusion protein of human ZFP64

Applications: ELISA, IHC

Recommended Dilutions: IHC: 25-100; ELISA: 5000-10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human, Mouse


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Epigenetics and Nuclear Signaling

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human esophagus cancer tissue using 218044(ZFP64 Antibody) at a dilution of 1/25(Nucleus).



Immunohistochemistry analysis of paraffin-embedded Human thyroid cancer tissue using 218044(Anti-ZFP64 Antibody) at a dilution of 1/25.