

ZNF202 RABBIT PAB

Cat.#: S218235

Product Name: Anti-ZNF202 Rabbit Polyclonal Antibody

Synonyms: ZSCAN42; ZKSCAN10

UNIPROT ID: O95125 (Gene Accession - BC013382)

Background: Enables DNA-binding transcription repressor activity, RNA polymerase II-specific and RNA polymerase II cis-regulatory region sequence-specific DNA binding activity. Involved in negative regulation of transcription by RNA polymerase II. Located in chromosome; nuclear body; and nucleolus. [provided by Alliance of Genome Resources, Apr 2022]

Immunogen: Fusion protein of human ZNF202

Applications: ELISA, IHC

Recommended Dilutions: IHC: 50-100; ELISA: 5000-10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Epigenetics and Nuclear Signaling

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human colorectal cancer tissue using 218235(ZNF202 Antibody) at a dilution of 1/75(Nucleus).