

ZNF559 RABBIT PAB

Cat.#: S212851

Product Name: Anti-ZNF559 Rabbit Polyclonal Antibody

Synonyms: NBLA00121

UNIPROT ID: Q9BR84 (Gene Accession - BC006436)

Background: May be involved in transcriptional regulation.

Immunogen: Fusion protein of human ZNF559

Applications: ELISA, IHC

Recommended Dilutions: IHC: 50-300; ELISA: 5000-10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Epigenetics and Nuclear Signaling

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin embedded Human esophagus cancer tissue using 212851(ZNF559 Antibody) at a dilution of 1/80(Nucleus).



In comparison with the IHC on the left, the same paraffin-embedded Human esophagus cancer tissue is first treated with the fusion protein and then with 212851(Anti-ZNF559 Antibody) at dilution 1/80.