

## HA TAG MAB

### Anti-HA Tag Mouse Monoclonal Antibody

**Cat. #:** 28003

**Size:** 100  $\mu$ L

**Symbol:** HA-tag

**Description:** Anti-HA Tag Mouse Monoclonal Antibody

**Background:** The HA-tag antibody recognizes any protein containing HA Tag.

**Immunogen:** HA

**Applications:** ELISA, WB, IF, IHC

**Recommended Dilutions:**

ELISA 1:2000-1:5000

WB 1:1000-1:5000

IF 1:100-1:200

**Concentration:** 1 mg/ml

**Host Species:** Mouse

**Format:** Liquid

**Clonality:** Monoclonal

**Isotype:** IgG


**Purity:** Purified from ascites

**Constituents:** PBS (without  $Mg^{2+}$  and  $Ca^{2+}$ ), pH 7.4, 150 mM NaCl, 50% glycerol

**Species Reactivity:** Recognizes protein contain HA.

**Storage Conditions:** Store at  $-20^{\circ}C$ . Avoid repeated freezing and thawing

**Western blot:**



Western blot: Recombinant HA tag fusion protein was diluted with HA mouse monoclonal antibody (1;1000), and the protein loading amount was 50ng, and Western blot analysis was performed.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---