

## CD45 MAB (APAMISTAMAB)

目录: 28148

产品名称: CD45 (Apamistamab) Monoclonal Antibody

基因符号: B220;CD45;CD45R;GP180;L-CA;LCA;LY5;T200

描述: CD45 (apamistamab ) Monoclonal Antibody

背景: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

经过测试的应用: Flow Cyt

推荐稀释度: Flow Cyt 1:100

种属反应性: Mus musculus

亚型: IgG1

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: Human CD45

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with 15 µg/mL CD45 (apamistamab biosimilar) mAb (28148) on Expi293 cells transfected with Human CD45 protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).

