

HUMAN APLNR FULL LENGTH PROTEIN

目录: 11113

产品名称: Human APLNR Full Length Protein

规格: 10 µg; 50 µg and 100 µg

基因符号: AGTRL1; APJ; APJR; HG11

Target: APLNR

UNIPROT ID: P35414

描述: Human APLNR Full Length Protein-Synthetic Nanodisc

背景: A member of the G protein-coupled receptor gene family. The encoded protein is related to the angiotensin receptor, but is actually an apelin receptor that inhibits adenylate cyclase activity and plays a counter-regulatory role against the pressure action of angiotensin II by exerting hypertensive effect. It functions in the cardiovascular and central nervous systems, in glucose metabolism, in embryonic and tumor angiogenesis and as a human immunodeficiency virus (HIV-1) coreceptor. Two transcript variants resulting from alternative splicing have been identified.

Species/Host: HEK293

Molecular Weight: The human full length APLNR protein has a MW of 42.7 kDa

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Protein Families: Druggable Genome, GPCR, Transmembrane

Protein Pathways: Neuroactive ligand-receptor interaction

ELISA assay to evaluate APLNR-Nanodisc 0.2µg Human APLNR-Nanodisc per well




Figure 1. Elisa plates were pre-coated with Flag Tag APLNR-Nanodisc (0.2µg/per well). Serial diluted Flag monoclonal antibody solutions were added, washed, and incubated with secondary antibody before Elisa reading. From above data, the EC50 for Flag monoclonal antibody binding with APLNR-Nanodisc is 2.278ng/ml.



Figure 2. Human APLNR-Nanodisc, Flag Tag on SDS-PAGE