

APOLIPOPROTEIN A1**Apolipoprotein A1 Monoclonal Antibody**

货号: 26025

基因符号: APOA1

描述: Apolipoprotein A1 小鼠单克隆抗体

背景: Apolipoprotein A1 is the major protein component of high-density lipoprotein (HDL) ("good cholesterol") in plasma. Apolipoprotein A1 helps to clear cholesterol from arteries.

免疫原: Recombinant full length protein of human Apolipoprotein A1

经过测试的应用: ELISA, WB, IHC, IP

推荐稀释度:

ELISA: 1:100-1:5000

WB: 1:100-1:2000

IHC: 1:50-1:100

浓度: 1 mg/ml

种属反应性: Mouse

形式: Liquid

克隆性: Monoclonal

亚型: IgG


纯化: Purified from ascites

Preservative: No

成分: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: Apolipoprotein A1 antibody recognizes Apolipoprotein A1 of primates.

储存条件: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing

Western blot:

WB: anti-Apo A1 mAb

Western blot analysis of Apolipoprotein A1 proteins.

