

## BIOTINYLATED ANTI-2B4 (DM69) RABBIT MAB

**Cat.#:** 28380

**Product Name:** Biotinylated Anti-2B4(DM69) Rabbit Monoclonal Antibody

**Synonyms:** CD244;2B4;SLAMF4;NKR2B4;NAIL;h2B4

**Description:** Biotinylated Anti-2B4 antibody(DM69) Rabbit Monoclonal Antibody

**Background:** This gene encodes a cell surface receptor expressed on natural killer (NK) cells (and some T cells) that mediate non-major histocompatibility complex (MHC) restricted killing. The interaction between NK-cell and target cells via this receptor is thought to modulate NK-cell cytolytic activity. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

**Applications:** ELISA; Flow Cyt

**Recommended Dilutions:** ELISA 1:5000-10000; Flow Cyt 1:100

**Host Species:** Rabbit

**Isotype:** Rabbit IgG

**Purification:** Purified from cell culture supernatant by affinity chromatography

**Species Reactivity:** Human 2B4

**Constituents:** Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

**Storage & Shipping:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).