

## **BIOTINYLATED ANTI-2B4 (DM70) RABBIT MAB**

目录: 28382

产品名称: Biotinylated 2B4(DM70) Rabbit Monoclonal Antibody

基因符号: CD244;2B4;SLAMF4;NKR2B4;NAIL;h2B4

描述: Biotinylated 2B4 antibody(DM70) Rabbit Monoclonal Antibody

背景: This gene encodes a cell surface receptor expressed on natural killer (NK) cells (and some T cells) that mediate non-major histocompatibility complex (MHC) restricted killing. The interaction between NK-cell and target cells via this receptor is thought to modulate NK-cell cytolytic activity. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

经过测试的应用: ELISA; Flow Cyt

推荐稀释度: ELISA 1:5000-10000; Flow Cyt 1:100

种属反应性: Rabbit

亚型: Rabbit IgG

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: Human 2B4

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).