

## **BIOTINYLATED ANTI-BCMA (DM4) RABBIT MAB**

目录: 28289

产品名称: Biotinylated BCMA(DM4) Rabbit Monoclonal Antibody

基因符号: TNFRSF17

描述: Biotinylated BCMA antibody(DM4) Rabbit Monoclonal Antibody

背景: The protein encoded by this gene is a member of the TNF-receptor superfamily. This receptor is preferentially expressed in mature B lymphocytes; and may be important for B cell development and autoimmune response. This receptor has been shown to specifically bind to the tumor necrosis factor (ligand) superfamily; member 13b (TNFSF13B:TALL-1:BAFF); and to lead to NF-kappaB and MAPK8:JNK activation. This receptor also binds to various TRAF family members; and thus may transduce signals for cell survival and proliferation. [provided by RefSeq; Jul 2008]

经过测试的应用: ELISA; Flow Cyt

推荐稀释度: Flow Cyt 1:100; IP 1:30

种属反应性: Rabbit

亚型: Rabbit IgG

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: Human BCMA

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).