

BIOTINYLATED ANTI-CD28 (DM63) RABBIT MAB

Cat.#: 28370

Product Name: Biotinylated Anti-CD28(DM63) Rabbit Monoclonal Antibody

Synonyms: CD28; Tp44

Description: Biotinylated Anti-CD28 antibody(DM63) Rabbit Monoclonal Antibody

Background: The protein encoded by this gene is essential for T-cell proliferation and survival; cytokine production; and T-helper type-2 development. Several alternatively spliced transcript variants encoding different isoforms have been found for this gene.

Applications: ELISA; Flow Cyt

Recommended Dilutions: ELISA 1:5000-10000; Flow Cyt 1:100

Host Species: Rabbit

Isotype: Rabbit IgG

Purification: Purified from cell culture supernatant by affinity chromatography

Species Reactivity: Human CD28

Constituents: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).