

BIOTINYLATED ANTI-TROP2 (DM75) RABBIT MAB

目录: 28392

产品名称: Biotinylated Trop2(DM75) Rabbit Monoclonal Antibody

基因符号: TACSTD2; GA733-1; MIS1; TROP2

描述: Biotinylated Trop2 antibody(DM75) Rabbit Monoclonal Antibody

背景: This intronless gene encodes a carcinoma-associated antigen. This antigen is a cell surface receptor that transduces calcium signals. Mutations of this gene have been associated with gelatinous drop-like corneal dystrophy.

经过测试的应用: ELISA; Flow Cyt

推荐稀释度: ELISA 1:5000-10000; Flow Cyt 1:100

种属反应性: Rabbit

亚型: Rabbit IgG

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: Human Trop2

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).