

CD19 (DM207) RABBIT MAB

目录: 28550

产品名称: CD19(DM207) Rabbit Monoclonal Antibody

基因符号: CD19;B4;CVID3;MGC12802

描述: CD19 antibody(DM207) Rabbit Monoclonal Antibody

背景: Lymphocytes proliferate and differentiate in response to various concentrations of different antigens. The ability of the B cell to respond in a specific; yet sensitive manner to the various antigens is achieved with the use of low-affinity antigen receptors. This gene encodes a cell surface molecule which assembles with the antigen receptor of B lymphocytes in order to decrease the threshold for antigen receptor-dependent stimulation.

经过测试的应用: ELISA; Flow Cyt

推荐稀释度: ELISA 1:5000-10000; Flow Cyt 1:100

种属反应性: Rabbit

亚型: Rabbit IgG

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: Human CD19

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with CD19 (DM207) on Expi293 cells transfected with human CD19 (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).