

CD40 LIGAND (DM164) RABBIT MAB

目录: 28507

产品名称: CD40 Ligand(DM164) Rabbit Monoclonal Antibody

基因符号: CD40LG;CD154;CD40L;HIGM1;IGM;IMD3;T-BAM;TNFSF5;TRAP;gp39

描述: CD40 Ligand antibody(DM164) Rabbit Monoclonal Antibody

背景: The protein encoded by this gene is expressed on the surface of T cells. It regulates B cell function by engaging CD40 on the B cell surface. A defect in this gene results in an inability to undergo immunoglobulin class switch and is associated with hyper-IgM syndrome.

经过测试的应用: ELISA; Flow Cyt

推荐稀释度: ELISA 1:5000-10000; Flow Cyt 1:100

种属反应性: Rabbit

亚型: Rabbit IgG

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: Human CD40 Ligand

成分: Lyophilized from sterile PBS, pH 7.4. 5 % - 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with CD40L (DM164) on Expi293 cells transfected with human CD40L (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram)

