

CD73 (DM204) RABBIT MAB

目录: 28547

产品名称: CD73(DM204) Rabbit Monoclonal Antibody

基因符号: CD73;NT5E;5'-Nucleotidase;5'-NT;NT5;NTE

描述: CD73 antibody(DM204) Rabbit Monoclonal Antibody

背景: The protein encoded by this gene is a plasma membrane protein that catalyzes the conversion of extracellular nucleotides to membrane-permeable nucleosides. The encoded protein is used as a determinant of lymphocyte differentiation. Defects in this gene can lead to the calcification of joints and arteries. Two transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq; Mar 2011]

经过测试的应用: ELISA; Flow Cyt

推荐稀释度: ELISA 1:5000-10000; Flow Cyt 1:100

种属反应性: Rabbit

亚型: Rabbit IgG

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: Human CD73

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with CD73 (DM204) on Expi293 cells transfected with human CD73 (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).