

50ML CENTRIFUGE TUBE FILTER, 0.22MM PES FILTER, STERILE

50mL Centrifuge Tube Filter Systems

- Eliminate a sample transfer step by filtering directly into a 50mL tube
- Polyethersulfone (PES) membrane for high flow and low binding
- One Centrifuge Tube Stand included with each case of 50mL Centrifuge Tube Filter Systems
- Centrifuge tube Stands are manufactured from polypropylene and are autoclavable, reusable and recyclable
- Centrifuge Tube Stand will hold one 50mL tube, three 15mL tubes and three microcentrifuge tubes (1.5mL)

Cat. #: CELLTREAT 227910

Pack Description: 1/Bag

Case Quantity: 12

Sterile: Yes

Non-pyrogenic: Yes

Primary Color: Clear

Secondary Color: Cap: Green

Primary Material: Tube: Polypropylene

Secondary Material: Filter: Polystyrene

Maximum Volume: 50mL

Working Volume: 50mL

Graduated: Yes

Graduation Increments (Subdivisions): 2.5mL

Max Force (RCF)g: 16,000

Pore Size: 0.22 μ m

RNase/DNase Free: Yes



Note: This product was over-stocked from CELLTREAT. It is 50% off CELLTREAT list price.