

CORNING® DISPOSABLE 100 MM ASEPTIC TRANSFER CAP FOR 5L PLASTIC ERLENMEYER FLASK, 0.2 MM VENT, MALE MPC, STERILE

Closed cell culture systems help reduce the risk of adventitious contamination during drug development and manufacturing. Corning offers closed systems accessories with vessels to help transfer liquid aseptically. The internal end is a rigid polypropylene dip tube designed to reach the bottom of the flask. The venting port ends in a 50mm diameter vent filter with a 0.2 µm pore PVDF membrane and a barbed fitting. The liquid transfer port uses a (1/4 OD) chemically resistant, heat sealable flexible tubing on the external side of the cap ending in a male plastic coupler (MPC) male disconnect with a sealing cap.

- Double bagged for GMP applications
- Animal component free vessels
- BSE/TSE compliant accessories
- SAL 10 \times sterilization
- Wide range of product offerings
- Customizable solutions

Product Number: Corning 11500

Qty./Pk: 1 / Pk

Qty./Cs: 4 / Cs

Brand: Corning®

Cap Style: Screw cap

Cap Material: Polypropylene

Cap Diameter: 100 mm

Cap Color: Orange

Color: Orange

Sterile: Yes



Note: This product was over-stocked from Corning. It is 70% off Corning list price.