

CORNING® 1L POLYCARBONATE ERLLENMEYER FLASK WITH DIP TUBE, WITH 0.2 MM PTFE FILTER, MLL

Corning Erlenmeyer Flasks are now available as closed system solutions. They are sterile and arrive individually double-bagged and preassembled. The preassembled transfer caps have two ports. One port ends in a 0.2 µm filter, and the other port is chemically resistant, heat sealable flexible tubing, ending in a male Luer lock. The dip tube reaches all the way to the bottom of the flask for easy aseptic transfer of your liquid handling processes.

- Made from optically clear polycarbonate
- Sterilized by gamma radiation and nonpyrogenic

Cat. #: Corning 11440

Qty./Pk: 1 / Pk

Qty./Cs: 4 / Cs

Brand: Corning®

Cap Material: Polypropylene

Cap Color: Orange

Sterile: Yes

Feature: Disposable

Capacity: 1000 mL

Flask Style: Erlenmeyer

Top Style: Threaded

Graduation Interval: 50 mL



Note: This product was over-purchased from Corning. It is 70% off Corning list price.