


**EGFR(T854A)****EGFR(T854A)****Cat. #:** 26481**Gene Symbol:** EGFR/HER1/ERBB1**Description:** Anti-EGFR(T854A) Mouse Monoclonal Antibody**Background:** EGFR is a receptor tyrosine kinase binding ligands of the EGF family and activating several signaling cascades to convert extracellular cues into appropriate cellular responses. It activates at least 4 major downstream signaling cascades including the RAS-RAF-MEK-ERK, PI3 kinase-AKT, PLCgamma-PKC and STATs modules. Defects in EGFR are associated with lung cancer (LNCR).**Immunogen:** A synthetic peptide from the internal region of EGFR which includes the mutation of T854A, human origin.**Applications:** ELISA, WB, IF, IHC**Recommended Dilutions:**

ELISA: 1:1000-1:5000

WB: 1:1000-1:5000

IHC: 1:50-1:200


**Concentration:** 1 mg/ml**Host Species:** Mouse**Format:** Liquid**Clonality:** Monoclonal**Isotype:****Purity:** Purified from ascites**Preservative:** No**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150 mM NaCl, 50% glycerol**Species Reactivity:** Recognizes EGFR (T854A), but not wild type EGFR from vertebrates.**Storage Conditions:** Store at -20°C. Avoid repeated freezing and thawing**Western blot:**



## Western blot analysis of recombinant EGFR (T854A) and wild type proteins.

Purified His-tagged EGFR (T854A) protein (lane 1) and corresponding wild type protein (amino acids 702-987, lane 2) were blotted with anti EGFR(T854A) mouse monoclonal antibody (Cat. #26481).

## Immunofluorescence:



## Immunofluorescence of cells expressing EGFR proteins with anti-EGFR (T854A) antibody.

HEK293T cells were transfected with pCDNA3-GFP-EGFR (WT) plasmid (left column) or pCDNA3-GFP EGFR(T854A) plasmid (right column), then fixed and stained with anti-EGFR (T854A) monoclonal antibody (Cat. #26481).



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---