

TWEAKR MAB (ENAVATUZUMAB)

目录: 28040

产品名称: TweakR (Enavatuzumab) Monoclonal Antibody

基因符号: TNFRSF12A;GF-inducible 14;N14;tweakR;D266

描述: TweakR (enavatuzumab) Monoclonal Antibody

背景: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

经过测试的应用: ELISA; Flow Cyt

推荐稀释度: ELISA 1:5000-10000; Flow Cyt 1:100

种属反应性: Humanized

亚型: IgG1

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: Human TweakR

成分: Lyophilized from sterile PBS, pH 7.4. 5% - 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).

Anti-TweakR (enavatuzumab biosimilar) mAb ELISA




Figure 1. ELISA plate pre-coated by 2 µg/ml (100 µl/well) Human TweakR, mFc tagged protein 11275 can bind TweakR Neutralizing antibody in a linear range of 0.64-80 ng/ml.




Figure 2. Flow cytometry analysis with 1 µg/mL TweakR (enavatuzumab) mAb (28040) on Expi293 cells transfected with Human TweakR protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).