

## CGRPR MAB (ERENUMAB)

**Cat.#:** 28157

**Product Name:** Anti-CGRPR(Erenumab) Monoclonal Antibody

**Synonyms:** CALCRL

**Description:** Anti-CGRPR(Erenumab ) Monoclonal Antibody

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Applications:** Flow Cyt

**Recommended Dilutions:** Flow Cyt 1:100

**Host Species:** Homo sapiens

**Isotype:** IgG2

**Purification:** Purified from cell culture supernatant by affinity chromatography

**Species Reactivity:** Human CGRPR

**Constituents:** Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

**Storage & Shipping:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with 1µg/mL Anti-CGRPR(Erenumab biosimilar) mAb (28157) on Expi293 cells transfected with Human CGRPR&RAMPI protein (Blue histogram) and Expi293 transfected with Human CGRPR protein (Red histogram) ) or Human RAMPI protein(Orange histogram).



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---