

## EYA2 PAB

### EYA2

货号: 21128

基因符号: EYA2

描述: EYA2 Rabbit Polyclonal Antibody

背景: EYA2 (eyes absent homolog 2) functions as a tyrosine phosphatase that specifically dephosphorylates Tyr-142 of histone H2AX. EYA2 promotes efficient DNA repair and coactivates SIX1. Together with SIX1 and DACH2, EYA2 seems to be involved in myogenesis.

免疫原: A synthetic peptide from the internal region of human EYA2.

经过测试的应用: ELISA, WB, IHC, IP

推荐稀释度:

ELISA: 1:1000-1:5000

WB: 1:100-1:1000

IHC: 1:100-1:1000

浓度: 0.1 mg/ml

种属反应性: Rabbit

形式:

克隆性: Polyclonal

亚型: IgG

纯化: Purified from serum


**Preservative:** No

成分: PBS (without  $Mg^{2+}$  and  $Ca^{2+}$ ), pH 7.4, 150 mM NaCl, 59% glycerol Species

种属反应性: EYA2 antibody recognizes EYA2 of vertebrates.

储存条件: Store at  $-20^{\circ}C$ . Avoid repeated freezing and thawing

## Immunoprecipitation/Western blot: EYA2 protein



### WB: Anti-EYA2 pAb

#### Western blot analysis of recombinant full length EYA2 protein.

Purified His-tagged EYA2 protein (货号10195) was blotted with EYA2 polyclonal antibody (Cat # 21128).