

FALCON® 48-WELL CLEAR FLAT BOTTOM NOT TREATED MULTIWELL CELL CULTURE PLATE, WITH LID, INDIVIDUALLY WRAPPED, STERILE, 50/CASE

Falcon® 48-well Clear Flat Bottom Not Treated Multiwell Cell Culture Plate, with Lid, Individually Wrapped, Sterile, 50/Case

Falcon® 48-well Clear Multiwell Plate is comprised of clear not treated polystyrene. Flat bottom, with low evaporation lid. Individually wrapped, sterile, 50/case. Unique labyrinth lid, condensation rings, and deep-well design control contamination, reduce evaporation, and minimize edge effects.

- Crystal-grade virgin polystyrene
- Nonpyrogenic
- Non-treated polystyrene
- Convenient, peel-open medical-style packaging for individual product presentation
- Individual trays are PET (Code 1) and recyclable
- Working volume of Media is 0.5–0.8 mL
- Working volume of Trypsin is 0.05–0.08 mL (Your lab protocol may call for another proteolytic enzyme.)

Cat. #: Falcon® 351178

Qty./Pk: 1 / Pk

Qty./Cs: 50 / Cs

Brand: Falcon®

Plate Format: 48-well

Well Shape: Round

Cell Growth Area: 0.75 cm²

Surface Treatment: Not Treated

Sterile: Yes

Lids Included: Yes



Note: This product was over-stocked from Corning. It is 50% off Corning list price.