

GA13 MAB

Gα₁₃ Monoclonal Antibody

货号: 26004

基因符号: Gna13

描述: Gα₁₃ 小鼠单克隆抗体

背景: Heterotrimeric G proteins are essential cellular signal transducers. Gα₁₃ is one of the G proteins that could mediate cell migration and angiogenesis. Other biochemical and physiological functions of Gα₁₃ are being explored.

免疫原: Recombinant full length Gα₁₃ protein

经过测试的应用: ELISA, WB, IHC

推荐稀释度: ELISA 1:1000–1:5000, WB 1:1000–1:2000

浓度: 1 mg/ml

种属反应性: Mouse

形式: Liquid

克隆性: Monoclonal

亚型: IgG

纯化: Purified from ascites


Preservative: No

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: Gα₁₃ monoclonal antibody recognizes Gα₁₃ of vertebrates.

储存条件: Store at -20°C. Avoid repeated freezing and thawing

Western blot:



Western blot analysis of Gα₁₃ proteins.

Gα₁₃, Gα₁₂ or Gα_q plasmid DNAs were transiently transfected into CHO cells (endogenous Gα₁₃ protein was down-regulated by Gα₁₃siRNA). Whole-cell lysates were used for western blots.