

ANTI GABA B RECEPTOR 1 NANOBODY**Anti GABA B Receptor 1 Nanobody**

货号: 26201

规格: 100 µL

基因符号: GABA B Receptor 1

描述: Anti GABA B Receptor 1 nanobody

背景: GABA B Receptor 1 is a subunit of the metabotropic γ-aminobutyric acid (GABA) B-type receptor. It is mainly responsible for binding and activating downstream G proteins and mediating signal transmission. GABA B Receptor 1 belongs to the C-type family of G protein-coupled receptors, regulates inhibitory neurotransmission, and is related to the onset of epilepsy, pain, addiction and other diseases. It is an important drug development target.

This product is a highly specific and high-affinity nanobody for GABA B Receptor 2. It recognizes the extracellular structural region of GABA B Receptor 1 and can be used in experiments such as GABA B Receptor 1 expression detection and immunoprecipitation. It is a good GABA B Receptor 1 detection tool.

免疫原: Recombinant monoclonal Nanobody, fused to Mouse Fc

经过测试的应用: Flow cytometry staining, immunostaining and Elisa, IP, Co-IP, IHC

Dilution: 1:400

Molecular Weight: 90kD

种属反应性: Mouse

形式: Liquid

克隆性: Monoclonal

Tag: His

纯化:: Purified from ascites

Preservative: No

成分: PBS, 40% Glycerol, 0.03% proclin300

Species: Human, Rat

储存条件:: Store at 4℃, 6 months; -20℃, 2 years

FACS GABA B Receptor 1 antibody : 

