

HTRA1 MAB (GALEGENIMAB)

目录: 28141

产品名称: HTRA1(Galegenimab) Monoclonal Antibody

基因符号: L56;Serine protease 11

描述: HTRA1(galegenimab) Monoclonal Antibody

背景: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

经过测试的应用: Flow Cyt

推荐稀释度: Flow Cyt 1:100

种属反应性: Humanized

亚型: IgG1

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: Human HTRA1

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with 15 µg/mL HTRA1(galegenimab biosimilar) mAb (28141) on Expi293 cells transfected with Human HTRA1 protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).

