

GA_o PAB

GA_o Polyclonal Antibody

货号: 21015

基因符号: Gnao

描述: GA_o Rabbit Polyclonal Antibody

背景: Heterotrimeric G proteins are critical cellular signal transducers. GA_o belongs to the Gi sub-family of G proteins that could mediate the inhibition of adenylyl cyclases. Other biochemical and physiological functions of GA_o proteins are being explored.

免疫原: Recombinant full length protein of GA_o.

经过测试的应用: ELISA, WB, IHC

推荐稀释度: ELISA 1:1000-1:5000

WB 1:100-1:2000

浓度: 1 mg/ml

种属反应性: Rabbit

形式: Liquid

克隆性: Polyclonal

亚型: IgG

纯化: Purified from serum


Preservative: No

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: GA_o antibody recognizes GA_o of vertebrates.

储存条件: Store at -20°C. Avoid repeated freezing and thawing

Western blot:



WB: anti-Gα_o pAb

Western blot analysis of Go proteins.

Whole-cell lysate of MEF cells or purified Gα_o protein were used for western blots with the polyclonal anti Gα_o antibody (Cat # 21015).