

GAQ PAB

Ga_q Polyclonal Antibody

货号: 21014

基因符号: Gnaq

描述: Gaq Rabbit Polyclonal Antibody

背景: Heterotrimeric G proteins are essential cellular signal transducers. Gaq is a member of the Gq subfamily of G proteins. Gaq is known to stimulate phospholipase C. Other biochemical and physiological functions of Gaq are being explored.

免疫原: Recombinant full length Gaq protein

经过测试的应用: ELISA, WB, IHC

推荐稀释度: ELISA 1:1000-1:5000

WB 1:100-1:2000

浓度: 1 mg/ml

种属反应性: Rabbit

形式: Liquid

克隆性: Polyclonal

亚型: IgG

纯化:: Purified from serum


Preservative: No

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: Gaq polyclonal antibody recognizes Gaq of vertebrates.

储存条件:: Store at -20°C. Avoid repeated freezing and thawing

Western blot:



WB: anti-Gα_q pAb

Western blot analysis of Gaq proteins.

MEF whole-cell lysates were used for western blot with Gaq polyclonal antibody (cat # 21014).



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
