

ANTI-GLP1R NANOBODY

Anti GLP1R Nanobody

货号: 26206

规格: 100 µL

基因符号: GLP1R

描述: GLP1R Nanobody

背景: GLP-1R (Glucagon-like peptide 1 receptor, glucagon-like peptide-1 receptor) is a class B, G-protein-linked seven-transmembrane domain protein found in the pancreas, lung, heart, kidney, Expressed in the hypothalamus and stomach. It is a receptor for glucagon-like peptide, a regulator of insulin secretion from pancreatic beta cells. This product is a nanobody with high specificity and affinity for GLP1R. It recognizes the extracellular structural region of GLP1R and can be used in experiments such as GLP1R expression detection and immunoprecipitation. It is a good GLP1R detection tool.

克隆性: Recombinant monoclonal Nanobody

经过测试的应用: Flow cytometry staining, immunostaining and Elisa, IP, Co-IP, IHC

Dilution: 1:3000

Molecular Weight: 14 kD

种属反应性: Human, Rat

形式: Liquid

成分: PBS, 40% Glycerol, 0.03% proclin300

种属反应性: Human, Mouse

储存条件: Store at 4℃, 6 months; -20℃, 2 years

FACS GLP1R antibody at 1/3000 dilution: 

Tag:

His,

Myc