

HUMAN GPR87 FULL LENGTH PROTEIN

目录: 11043

产品名称: Human GPR87 Full Length Protein

规格: 10 µg, 50 µg and 100 µg

基因符号: FKSG78; GPR95; KPG_002

Target: GPR87

UNIPROT ID: Q9BY21

描述: Human GPR87 full length protein-synthetic nanodisc

背景: This gene encodes a G protein-coupled receptor and is located in a cluster of G protein-couple receptor genes on chromosome 3. The encoded protein has been shown to be overexpressed in lung squamous cell carcinoma (PMID:18057535) and regulated by p53 (PMID:19602589). [provided by RefSeq, Nov 2011]

Species/Host: HEK293

Molecular Weight: The human full length GPR87 protein has a MW of 41.3 kDa

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.