

## HIS TAG (DM1) RABBIT MAB

目录: 28287

产品名称: His Tag(DM1) Rabbit Monoclonal Antibody

基因符号: his tag

描述: His Tag Antibody(DM1) Rabbit Monoclonal Antibody

背景: The rabbit immune system generates antibody diversity and optimizes affinity. His Tag Antibody(DM1); Rabbit mAb is specific to six histidine tags placed at C-terminal and internal regions of fusion proteins. This antibody can greatly improve the effectiveness of several different kinds of immunoassays; helping researchers identify; detect; and purify polyhistidine fusion proteins in bacteria; insect cells; and mammalian cells.

经过测试的应用: ELISA; WB

推荐稀释度: WB 1:1000

种属反应性: Rabbit

亚型: Rabbit IgG

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: Detects proteins containing accessible consecutive histidine regions. The antibody detects His tags localized at the amino- or carboxyl- terminus. His Tag

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Western blot analysis of his-tag fusion protein (70kDa, 20ng) using His-Tag antibody (28287) at 1:1000 dilution. Secondary antibody: HRP Goat Rabbit IgG (H+L) at 1:5000 dilution.