

## HUMAN 2B4 PROTEIN, HFC TAG

**Cat.#:** 11287

**Product Name:** Human 2B4 Protein

**Size:** 10 µg, 50 µg and 100 µg

**Synonyms:** CD244;2B4;SLAMF4;NKR2B4;NAIL;h2B4

**Target:** 2B4

**UNIPROT ID:** Q9BZW8

**Description:** Recombinant human 2B4 protein with C-terminal human Fc tag

**Background:** This gene encodes a cell surface receptor expressed on natural killer (NK) cells (and some T cells) that mediate non-major histocompatibility complex (MHC) restricted killing. The interaction between NK-cell and target cells via this receptor is thought to modulate NK-cell cytolytic activity. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

**Species/Host:** HEK293

**Molecular Weight:** The protein has a predicted molecular mass of 47.7 kDa after removal of the signal peptide.

**Molecular Characterization:** 2B4(Cys22-Ala221) hFc(Glu99-Ala330)

**Purity:** The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

**Storage & Shipping:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human 2B4, hFc Tag on SDS-PAGE under reducing condition.

### Human 2B4, hFc Tagged protein ELISA

0.2 µg of 2B4, hFc Tagged protein per well




Figure 2. ELISA plate pre-coated by 2 µg/ml (100 µl/well) Human CD48, mFc-His tagged protein 11144 can bind Human 2B4, hFc tagged protein (11287) in a linear range of 62.5-4000 ng/ml.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---