

## HUMAN ABCB5 PROTEIN, HFC TAG

**Cat.#:** 11660

**Product Name:** Human ABCB5 Protein

**Size:** 10 µg, 50 µg and 100 µg

**Synonyms:** ABCB5alpha;ABCB5beta;EST422562

**Target:** ABCB5

**UNIPROT ID:** Q2M3G0

**Description:** Recombinant human ABCB5 protein with C-terminal human Fc tag

**Background:** ABCB5 belongs to the ATP-binding cassette (ABC) transporter superfamily of integral membrane proteins. These proteins participate in ATP-dependent transmembrane transport of structurally diverse molecules ranging from small ions, sugars, and peptides to more complex organic molecules (Chen et al., 2005 [PubMed 15760339]).  
[supplied by OMIM, Mar 2008]

**Species/Host:** HEK293

**Molecular Weight:** The protein has a predicted molecular mass of 43.34 kDa after removal of the signal peptide.

**Molecular Characterization:** ABCB5(Tyr757-Val827) hFc(Glu99-Ala331)

**Purity:** The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

**Storage & Shipping:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.