

**HUMAN ADAM28 PROTEIN, HIS TAG****Cat.#:** 11511**Product Name:** Human ADAM28 Protein**Size:** 10 µg, 50 µg and 100 µg**Synonyms:** eMDC II;MDC-L**Target:** ADAM28**UNIPROT ID:** Q9UKQ2**Description:** Recombinant human ADAM28 protein with C-terminal 6xHis tag**Background:** This gene encodes a member of the ADAM (a disintegrin and metalloprotease domain) family. Members of this family are membrane-anchored proteins structurally related to snake venom disintegrins, and have been implicated in a variety of biological processes involving cell-cell and cell-matrix interactions, including fertilization, muscle development, and neurogenesis. The protein encoded by this gene is a lymphocyte-expressed ADAM protein. This gene is present in a gene cluster with other members of the ADAM family on chromosome 8. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2015]**Species/Host:** HEK293**Molecular Weight:** The protein has a predicted molecular mass of 73.7 kDa after removal of the signal peptide. The apparent molecular mass of ADAM28-His is approximately 55-70 kDa and 70-100 kDa due to glycosylation.**Molecular Characterization:** ADAM28(Ile19-His664) 6xHis tag**Purity:** The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.**Storage & Shipping:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human ADAM28 Protein, His Tag on SDS-PAGE under reducing condition.