

HUMAN CD1A PROTEIN, HFC TAG

Cat.#: 11635

Product Name: Human CD1A Protein; HFc Tag

Size : 10 µg; 50 µg and 100 µg

Synonyms: R4; T6; CD1; FCB6; HTA1

Target: CD1A

UNIPROT ID: P06126

Description: Recombinant Human CD1A Protein With C-Terminal Human Fc Tag

Background: This gene encodes a member of the CD1 family of transmembrane glycoproteins, which are structurally related to the major histocompatibility complex (MHC) proteins and form heterodimers with beta-2-microglobulin. The CD1 proteins mediate the presentation of primarily lipid and glycolipid antigens of self or microbial origin to T cells. The human genome contains five CD1 family genes organized in a cluster on chromosome 1. The CD1 family members are thought to differ in their cellular localization and specificity for particular lipid ligands. The protein encoded by this gene localizes to the plasma membrane and to recycling vesicles of the early endocytic system. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2016]

Species/Host: HEK293

Molecular Weight: The protein has a predicted molecular mass of 58.1 kDa after removal of the signal peptide.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human CD1A Protein, hFc Tag on SDS-PAGE under reducing condition.