

HUMAN CD26 PROTEIN, HFC TAG**Cat.#:** 11662**Product Name:** Human CD26 Protein**Size:** 10 µg, 50 µg and 100 µg**Synonyms:** ADABP;ADCP2;CD26;DPPIV;TP103**Target:** ADABP**UNIPROT ID:** P27487**Description:** Recombinant human CD26 protein with C-terminal human Fc tag**Background:** The DPP4 gene encodes dipeptidyl peptidase 4, which is identical to adenosine deaminase complexing protein-2, and to the T-cell activation antigen CD26. It is an intrinsic type II transmembrane glycoprotein and a serine exopeptidase that cleaves X-proline dipeptides from the N-terminus of polypeptides. Dipeptidyl peptidase 4 is highly involved in glucose and insulin metabolism, as well as in immune regulation. This protein was shown to be a functional receptor for Middle East respiratory syndrome coronavirus (MERS-CoV), and protein modeling suggests that it may play a similar role with SARS-CoV-2, the virus responsible for COVID-19. [provided by RefSeq, Apr 2020]**Species/Host:** HEK293**Molecular Weight:** The protein has a predicted molecular mass of 116.71 kDa after removal of the signal peptide.**Molecular Characterization:** ADABP(Asn29-Pro766) hFc(Glu99-Ala330)**Purity:** The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.**Storage & Shipping:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.