

## HUMAN CD72 PROTEIN, HFC TAG

**Cat.#:** 11555

**Product Name:** Human CD72 Protein

**Size:** 10 µg, 50 µg and 100 µg

**Synonyms:** Lyb-2

**Target:** CD72

**UNIPROT ID:** P21854

**Description:** Recombinant human CD72 protein with N-terminal human Fc tag

**Background:** Plays a role in B-cell proliferation and differentiation.[UniProtKB/Swiss-Prot Function]

**Species/Host:** HEK293

**Molecular Weight:** The protein has a predicted molecular mass of 54.2 kDa after removal of the signal peptide. The apparent molecular mass of hFc-CD72 is approximately 55-70 kDa due to glycosylation.

**Molecular Characterization:** hFc(Glu99-Ala330) CD72(Arg117-Asp359)

**Purity:** The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

**Storage & Shipping:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human CD72 Protein, hFc Tag on SDS-PAGE under reducing condition.