

HUMAN CDH3 PROTEIN, HIS TAG

Cat.#: 11588

Product Name: Human CDH3 Protein

Size: 10 µg, 50 µg and 100 µg

Synonyms: CDHP;HJMD;PCAD

Target: CDH3

UNIPROT ID: P22223

Description: Recombinant human CDH3 protein with C-terminal 6xHis tag

Background: This gene encodes a classical cadherin of the cadherin superfamily. Alternative splicing results in multiple transcript variants, at least one of which encodes a preproprotein that is proteolytically processed to generate the mature glycoprotein. This calcium-dependent cell-cell adhesion protein is comprised of five extracellular cadherin repeats, a transmembrane region and a highly conserved cytoplasmic tail. This gene is located in a gene cluster in a region on the long arm of chromosome 16 that is involved in loss of heterozygosity events in breast and prostate cancer. In addition, aberrant expression of this protein is observed in cervical adenocarcinomas. Mutations in this gene are associated with hypotrichosis with juvenile macular dystrophy and ectodermal dysplasia, ectrodactyly, and macular dystrophy syndrome (EEMS). [provided by RefSeq, Nov 2015]

Species/Host: HEK293

Molecular Weight: The protein has a predicted molecular mass of 60.9 kDa after removal of the signal peptide. The apparent molecular mass of CDH3-His is approximately 70–100 kDa due to glycosylation.

Molecular Characterization: CDH3(Ser18-Gly654) 6×His tag

Purity: The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human CDH3 Protein, His Tag on SDS-PAGE under reducing condition.