

HUMAN CEACAM8 PROTEIN, HIS TAG

目录: 11814

产品名称: Human CEACAM8 Protein

规格: 10 µg, 50 µg and 100 µg

基因符号: CD66b;CD67;CGM6;NCA-95

Target: CEACAM8

UNIPROT ID: P31997

描述: Recombinant human CEACAM8 protein with C-terminal 6xHis tag

背景: Cell surface glycoprotein that plays a role in cell adhesion in a calcium-independent manner (PubMed:8776764, PubMed:2022629, PubMed:11590190). Mediates heterophilic cell adhesion with other carcinoembryonic antigen-related cell adhesion molecules, such as CEACAM6 (PubMed:8776764, PubMed:2022629, PubMed:11590190). Heterophilic interaction with CEACAM8 occurs in activated neutrophils (PubMed:8776764).[UniProtKB/Swiss-Prot Function]

Species/Host: HEK293

Molecular Weight: The protein has a predicted molecular mass of 32.2 kDa after removal of the signal peptide. The apparent molecular mass of CEACAM8-His is approximately 35–70 kDa due to glycosylation.

Molecular Characterization: CEACAM8(Gln35–Ser319) 6×His tag

纯化: The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at –20°C to –80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at –80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human CEACAM8 Protein, His Tag on SDS-PAGE under reducing condition.