

HUMAN FGF-8F (N-6HIS) PROTEIN

目录: 12086

产品名称: Human FGF-8F (N-6His) Protein

规格: 10 µg, 50 µg and 100 µg

基因符号: Fibroblast growth factor 8; Androgen-induced growth factor; Heparin-binding growth factor 8; AIGF; HBGF-8; FGF-8B

Target: FGF-8f

UNIPROT ID: P55075-4

描述: Recombinant Human Fibroblast Growth Factor 8f is produced by our E.coli expression system and the target gene encoding Gln23-Arg244 is expressed with a 6His tag at the N-terminus.

背景: Fibroblast growth factor 8 (FGF-8) is a member of the fibroblast growth factor family. It is discovered as a growth factor essential for the androgen-dependent growth of mouse mammary carcinoma cells. Mouse FGF-8b shares 100% aa identity with human FGF-8b. FGF-8 is widely expressed during embryogenesis, and mediates epithelial-mesenchymal transitions. It plays an important role in the regulation of embryonic development, cell proliferation, cell differentiation and cell migration. It is required for normal brain, eye, ear, limb development during embryogenesis and normal development of the gonadotropin-releasing hormone (GnRH) neuronal system.

Species/Host: E.coli

Molecular Weight: 27.7 KDa

Molecular Characterization: Not available

纯化: Greater than 95% as determined by reducing SDS-PAGE.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Greater than 95% as determined by reducing SDS-PAGE.