

HUMAN FMNL1 PROTEIN, HFC TAG

目录: 11757

产品名称: Human FMNL1 Protein

规格: 10 µg, 50 µg and 100 µg

基因符号: C17orf1;C17orf1B;FHOD4;FMNL;KW-13

Target: FMNL1

UNIPROT ID: O95466

描述: Recombinant human FMNL1 protein with C-terminal human Fc tag

背景: This gene encodes a formin-related protein. Formin-related proteins have been implicated in morphogenesis, cytokinesis, and cell polarity. An alternative splice variant has been described but its full length sequence has not been determined. [provided by RefSeq, Jul 2008]

Species/Host: HEK293

Molecular Weight: The protein has a predicted molecular mass of 146.41 kDa after removal of the signal peptide.

Molecular Characterization: FMNL1 (His2-Leu1100) hFc(Glu99-ALA330)

纯化: The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.