

## HUMAN GITR PROTEIN, MFC-HIS TAG

**Cat.#:** 11147

**Product Name:** Human GITR Protein

**Size:** 10 µg, 50 µg and 100 µg

**Synonyms:** AITR;GITR;TNFRSF18;CD357

**Target:** GITR

**UNIPROT ID:** Q9Y5U5

**Description:** Recombinant human GITR protein with C-terminal mouse Fc and 6xHis tag

**Background:** This gene encodes a member of the TNF-receptor superfamily. The encoded receptor has been shown to have increased expression upon T-cell activation, and it is thought to play a key role in dominant immunological self-tolerance maintained by CD25(+)CD4(+) regulatory T cells. Knockout studies in mice also suggest the role of this receptor is in the regulation of CD3-driven T-cell activation and programmed cell death. Three alternatively spliced transcript variants of this gene encoding distinct isoforms have been reported.

**Species/Host:** HEK293

**Molecular Weight:** The protein has a predicted molecular mass of 50–52 kDa after removal of the signal peptide.

**Molecular Characterization:** GITR(Gln26-Pro162) mFc(Pro99-Lys330) 6×His

**Purity:** The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

**Storage & Shipping:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human GITR Protein, mFc-His Tag on SDS-PAGE under reducing condition.




Figure 2. HEK293 cell line transfected with irrelevant protein (B) and human GITR Ligand (A) were surface stained with Human GITR, mFc-His tagged protein (11147) 1µg/ml followed by Alexa 488-conjugated anti-mouse IgG secondary antibody.