

HUMAN IFN ALPHA4 (C-6HIS) PROTEIN

Cat.#: 12070

Product Name: Human IFN Alpha4 (C-6His) Protein

Size: 10 µg, 50 µg and 100 µg

Synonyms: Interferon alpha-4; Interferon alpha-4B; Interferon alpha-76; Interferon alpha-M1; IFNA4

Target: IFN

UNIPROT ID: P05013

Description: Recombinant Human Interferon Alpha-4 is produced by our Mammalian expression system and the target gene encoding Cys24-Asp189 is expressed with a 6His tag at the C-terminus.

Background: Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase. They allow for communication between cells to trigger the protective defenses of the immune system that eradicate pathogens or tumors. They are typically divided among three IFN classes: Type I, Type II and Type III. IFNA4 is a secreted protein and produced by macrophages. Two variants of IFNA4 (IFNA4a and IFNA4b) are known, which differ from each other by changes in their coding regions at nucleotide positions 220 and 410.

Species/Host: HEK293

Molecular Weight: 20.4 KDa

Molecular Characterization: Not available

Purity: Greater than 95% as determined by reducing SDS-PAGE.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Greater than 95% as determined by reducing SDS-PAGE.