

## HUMAN IFNAR1 PROTEIN, HIS TAG

**Cat.#:** 11692

**Product Name:** Human IFNAR1 Protein

**Size:** 10 µg, 50 µg and 100 µg

**Synonyms:** AVP;IFN- $\alpha$ -REC;IFNAR;IFNBR;IFRC

**Target:** IFNAR1

**UNIPROT ID:** P17181

**Description:** Recombinant Human IFNAR1 Protein with C-terminal 6xHis tag

**Background:** The protein encoded by this gene is a type I membrane protein that forms one of the two chains of a receptor for interferons alpha and beta. Binding and activation of the receptor stimulates Janus protein kinases, which in turn phosphorylate several proteins, including STAT1 and STAT2. The protein belongs to the type II cytokine receptor family and functions as an antiviral factor. [provided by RefSeq, Jul 2020]

**Species/Host:** HEK293

**Molecular Weight:** The protein has a predicted molecular mass of 48.0 kDa after removal of the signal peptide. The apparent molecular mass of IFNAR1-His is approximately 70-100 kDa due to glycosylation.

**Molecular Characterization:** IFNAR1(Lys28-Lys436) 6×His Tag

**Purity:** The purity of the protein is greater than 90% as determined by SDS-PAGE and Coomassie blue staining.

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

**Storage & Shipping:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human IFNAR1 Protein, His Tag on SDS-PAGE under reducing condition.

**Human IFNAR1, His Tagged protein ELISA**  
0.1  $\mu$ g of Human IFNAR1, His tagged protein per well




Figure 2. ELISA plate pre-coated by 1  $\mu$ g/ml (100  $\mu$ l/well) Human IFNAR1 Protein, His Tag (11692) can bind Anti-IFNAR1 Neutralizing antibody 28115 in a linear range of 0.61-39.06 ng/mL.