

## HUMAN IL23A AND IL12B HETERODIMER PROTEIN, HFC TAG AND HIS TAG

**Cat.#:** 11896

**Product Name:** Human IL23A And IL12B Heterodimer Protein; HFc Tag And His Tag

**Size:** 10 µg; 50 µg and 100 µg

**Synonyms:** IL-23;IL-23A;IL23P19;P19;SGRF and CLMF;CLMF2;IL-12B;IMD28;IMD29;NKSF;NKSF2

**Target:** IL23A andamp;IL12B

**UNIPROT ID:** Q9NPF7;P29460

**Description:** Recombinant Human IL23A Protein With C-Terminal Human Fc Tag And Human IL12B Protein With C-Terminal 6×His Tag

**Background:** Interleukin-23 subunit alpha (IL-23 alpha) can associates with IL12B to form the IL-23 interleukin, a heterodimeric cytokine which functions in innate and adaptive immunity. IL-23 may constitute with IL-17 an acute response to infection in peripheral tissues. IL-23 binds to a heterodimeric receptor complex composed of IL12RB1 and IL23R, activates the Jak-Stat signaling cascade, stimulates memory rather than naive T-cells and promotes production of proinflammatory cytokines. IL-23 induces autoimmune inflammation and thus may be responsible for autoimmune inflammatory diseases and may be important for tumorigenesis.

**Species/Host:** HEK293

**Molecular Weight:** The protein has a predicted molecular mass of 44.8 and 35.5 kDa after removal of the signal peptide. The apparent molecular mass of IL23A-hFc and IL12B-His is approximately 35-55 kDa due to glycosylation.

**Molecular Characterization:** IL23A(Arg20-Pro189) hFc(Glu99-Ala330)and IL12B(Ile23-Ser328) 6×His tag

**Purity:** The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HC, 150 mM NaC, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

**Storage & Shipping:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human IL23A & IL12B Heterodimer Protein, hFc Tag & His Tag on SDS-PAGE under reducing condition.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---