

HUMAN ITGAV-HIS AND ITGB1-HFC HETERODIMER PROTEIN

Cat.#: 12246

Product Name: Human ITGAV-His And ITGB1-HFc Heterodimer Protein

Size: 10 µg, 50 µg and 100 µg

Synonyms: Integrin alpha V beta 1;ITGAVandITGB1

Target: ITGAV andamp;ITGB1

UNIPROT ID: ITGAV;ITGB1

Description: Recombinant heterodimer protein containing both human ITGAV Protein with C-terminal 6xHis tag and human ITGB1 protein with C-terminal human Fc tag

Background: Integrin alpha-5/beta-1 is a receptor for fibrinogen. Integrin alpha-1/beta-1, alpha-2/beta-1, alpha-6/beta-1 and alpha-7/beta-1 are receptors for laminin. Integrin alpha-4/beta-1 is a receptor for VCAM1. It recognizes the sequence Q-I-D-S in VCAM1. Integrin alpha-9/beta-1 is a receptor for VCAM1, cytotactin and osteopontin. It recognizes the sequence A-E-I-D-G-I-E-L in cytotactin. Integrin alpha-V/beta-1 is also a receptor for vitronectin. Beta-1 integrins recognize the sequence R-G-D in a wide array of ligands. Isoform 2 interferes with isoform 1 resulting in a dominant negative effect on cell adhesion and migration (in vitro). When associated with alpha-7/beta-1 integrin, regulates cell adhesion and laminin matrix deposition.

Species/Host: HEK293

Molecular Weight: The protein has a predicted molecular mass of 107.1 and 104.5 kDa after removal of the signal peptide. The apparent molecular mass of ITGAV-His and ITGB1-hFc is approximately 130-250 kDa due to glycosylation.

Molecular Characterization: ITGAV(Phe31-Val992) 6xHis tag - ITGB1(Gln21-Asp728) hFc(Glu99-Ala330)

Purity: The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human ITGAV & ITGB1 Heterodimer Protein, His Tag & hFc Tag on SDS-PAGE under reducing condition.