

HUMAN LAMP PROTEIN, HFC TAG

Cat.#: 11700

Product Name: Human LAMP Protein

Size: 10 µg, 50 µg and 100 µg

Synonyms: CD208;DC-LAMP;DC LAMP;DCLAMP;LAMP;LAMP-3;TSC403

Target: LAMP

UNIPROT ID: Q13449

Description: Recombinant human LAMP protein with C-terminal human Fc tag

Background: Dendritic cells (DCs) are the most potent antigen-presenting cells. Immature DCs efficiently capture antigens and differentiate into interdigitating dendritic cells (IDCs) in lymphoid tissues that induce primary T-cell responses (summary by de Saint-Vis et al., 1998 [PubMed 9768752]). [supplied by OMIM, Dec 2010]

Species/Host: HEK293

Molecular Weight: The protein has a predicted molecular mass of 45.9 kDa after removal of the signal peptide.

Molecular Characterization: LAMP(Val29-Asn315) hFc(Glu99-Ala330)

Purity: The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.